


## ÜBER | ABOUT CERAMIX CUT © ★

KERAMIKBESCHICHTETE SCHERENEINSÄTZE | CERAMIC COATED SCISSOR INSERTS

### FAKTEN & VORTEILE | FACTS & FEATURES



Innovative Keramikbeschichtung | Innovative ceramic coating



- ▶ verschleißfest
  - ▶ verminderter Abrieb
  - ▶ kosteneffizient durch langzeitstabile Oberfläche
  - ▶ minimale optische Reflexion
  - ▶ robust
  - ▶ rost- und korrodierungsfrei
  - ▶ temperatur- und kratzbeständig
  - ▶ hervorragende Schneidfähigkeit
  - ▶ 3 bis 5fache Härte gegenüber gehärtetem Metall
  - ▶ glatte, dichte Struktur - leicht zu reinigende Oberfläche
  - ▶ beständiger gegenüber Reinigungs- und Desinfektionsmittel
  - ▶ verbesserte Gleiteigenschaft, dadurch weicher und gleichmäßiger Gang der Scherenblätter
  - ▶ geeignet zur Verwendung mit XPress Lock ©
- ▶ wear-resistant
  - ▶ reduced abrasion
  - ▶ cost-effective due to long-term stable surface
  - ▶ minimal optical reflection
  - ▶ solid
  - ▶ stainless and non-corrosive
  - ▶ temperature- and scratch-resistant
  - ▶ outstanding cutting ability
  - ▶ 3 to 5 times harder than hardened metal
  - ▶ smooth and dense surface - easily cleanable
  - ▶ more constant to detergents and disinfectants
  - ▶ improved sliding characteristics, therefore slighter movement of the scissor blades
  - ▶ suitable for the use with XPress Lock ©